

Drytop dryer cabinets

For all your flexible scopes



Your specialist in ENT treatment units, custom-made installations, patient chairs, dryer cabinets and digital storage systems. With more than 20 years of experience we guarantee top service, high quality and an ergonomic workplace. MediTop builds and develops all its products entirely in-house.



Drytop dryer cabinets

- Long endoscopes hang freely and off the ground
- Large automatic glass doors For good visibility , handsfree
- Own storage system for monitoring the storage time , pressure and flow
- Visual and audible alarm with irregularities
- Hepa filtered overpressure
- Connection with hygiene tracker optional
- Connection with current hospital systems optional
- accessible to all makes flexible scopes , whether or not with biopsy channel
- 8 to 20 scopes per module(depending the discipline)
- 2 scope modules connectable per technical part

we are pleased to provide additional information.

Dimensions:

